

SCHILLER
A M E R I C A S



TMX-428

High-quality, low-maintenance treadmill

High-quality, low-maintenance

SCHILLER

A M E R I C A S

MEDICAL TREADMILL

This heavy-duty, medical treadmill is specifically designed to be interfaced and controlled by Schiller Stress Test and a wide variety of Cardiac and Pulmonary Stress Testing Systems.

The TMX428 is distinguished by its smooth, quiet and reliable, while centering the design on the absolute safety and comfort of the patient. The treadmill is self-calibrating, virtually maintenance-free and the safety features include a low profile deck, .1MPH minimum operating speed, incredibly smooth belt operation and multiple emergency stop switches.

FEATURES

- Smooth, quiet, zero-start, 4HP (6HP peak) hi-torque drive
- 500 lb. patient capacity
- Smooth, quiet linear actuator elevation system to 25%
- Low-to-the ground, cushioned deck for patient comfort and safety
- 22" x 63" (56cm x 160cm) running belt and safety handrails
- Belt alignment system eliminates off-center running belts
- Exceptionally accurate, self-calibrating speed and elevation
- Prominent emergency stop button and patient safety tether
- Minimum maintenance



TMX-428

Technical Specifications

SCHILLER
A M E R I C A S

- 1.- Running Surface
- 2.- Adult/Pediatric Handrails
- 3.- Emergency Stop Button
- 4.- Emergency Stop Cord
- 5.- Elevation Support



TMX-428

Technical Specifications

SCHILLER
A M E R I C A S



GENERAL SPECIFICATIONS

Base dimension	84" x 36" (L x W)
Weight	480 lbs.
Max. Patient weight	500 lbs. / 227 kg
Safe, Low-profile deck	7" (18 cm) from floor
Running surface	22" x 63" (56 x 160 cm)
Self-aligning running belt	Yes
Cushioned and reversible running deck	Yes
All steel chassis	Yes
Powder Coat finish	Yes
Emergency stop and patient tether	Yes
Space Requirements	33"x 78.5" (84 x 199 cm)
Stress Test Interface	RS-232

Drive System

Heavy-duty, 4 HP (6 HP peak) drive

Power supply 208-240V, 50/60 Hz, dedicated 15A service
110V, 50/60 Hz, dedicated 20A service

Speed Range

Self-calibrating Yes

Elevation Range 0-25%, Self-calibrating

Other Options

Emergency stop button left or right side

Adjustable front handrail



TMX-428

SCHILLER
A M E R I C A S

Call us: 1-888-845-8455
Visit us at schillerservice.com
10903 NW 33rd St, Doral, FL 33172, USA